



**MEASUREMENT PARAMETERS**



FRIABILITY



ABRASION

**ADVANTAGES**

- + Compact stainless steel housing
- + Quick drum assembly
- + Automatic emptying

**OPTIONS**

- Friability test drum
- Abrasion test drum

**STANDARD DEVICE**

# AE 1 + 2

## Friability and abrasion tester

The abrasion and friability testing device consists of one or two test stations. If required, test drums can be equipped to test for wear due to rolling and falling or to test the abrasion of tablets, tablet cores, granules etc.

The tester meets all the latest requirements of the European & US Pharmacopoeia, GMP and GLP.

- Quick drum assembly
- User-friendly membrane keypad
- GMP-compliant stainless steel housing

MODEL	L x W x H (mm)	Weight (kg)
<b>AE 1</b>	320 x 300 x 250	6.5
<b>AE 2</b>	320 x 430 x 250	7



Friability and abrasion testing device AE 2 with 2 test drums