



H - SERIES EMMI-H30



The unit is supplied with lid and basket.



TECHNICAL DATA **EMMI-H30**

Description:

The Emmi-H30 offers plenty of space and a convincing cleaning result. Due to the exact shape of the basket, the lid is ideally placed on the device and kept heated for longer.

Dimensions:

External dimensions (L x W x H):
320 x 180 x 240 mm
Inside of Bath (L x W x H): 238 x 138 x 100 mm
Inside of basket (L x W x H): 205 x 105 x 70 mm

Content:

Tank Capacity: 3.0 Liters
Working Capacity: 2.4 Liters

Equipment:

Vibrator pan: stainless steel 1.4301
Housing: stainless steel 1.4301
Degree of protection: IP 32
Drain valve: Yes
Timer: 1-60 min / continuous operation
Heating: 20-80 °C

Technical data:

Oscillating system:
3PZT wide beam vibratory system
HF Generator: Overload protected, constant
Power operating mode: simple
Half-wave operating frequency: approx. 40 kHz
Ultrasonic power: max. 180 watts
Total wattage: 360 W
Heating power: 180 watts
Line voltage: 50 / 60 Hz
Power connection: 220 - 240 V

Accessories:

Stainless steel lid
«noise-resistant»
Stainless steel basket
«protects cleaning material and vibrating pan»