



## D - SERIES EMMI-D60



The unit is supplied with lid and basket.



### TECHNICAL DATA **EMMI-D60**

**Description:**

The Emmi-D60 is an ultrasonic device for the medium and higher Cleaning segment. Due to the large supply of space and increased cleaning performance, the D60 offers ideal cleaning results.

**Dimensions:**

External dimensions (L x W x H):  
380 x 175 x 240 mm  
Inside of Bath (L x W x H): 302 x 152 x 150 mm  
Inside of basket (L x w x H): 250 x 110 x 110 mm

**Capacity:**

Tank Capacity: 6.0 Liters  
Working Capacity: 4.7 Liters

**Equipment:**

Vibrator pan: stainless steel 1.4301  
Housing: stainless steel 1.4301  
Degree of protection: IP 32  
Drain valve: Yes  
Timer: 1-60 min / continuous operation  
Heating: 20-80 °C  
Ultrasonic power control: 50 / 100%

Sweep technique: yes (increase of the normal ultrasonic power in case of particularly strong soiling)

Degas function: yes (degassing of the tap water (e.g. oxygen) before starting results in an optimal ultrasonic power)

**Technical data:**

Oscillating system:  
4PZT wide beam vibratory system  
HF Generator: Overload protected, constant  
Power operating mode: simple  
Half-wave operating frequency: approx. 37 kHz  
Ultrasonic power: max. 320 watts  
Total power in Watt: 570 W  
Heating power: 250 Watt  
Line voltage: 50 / 60 Hz  
Power connection: 220 - 240 V

**Accessories:**

Stainless steel lid «noise-resistant»  
Stainless steel basket  
«protects cleaning material and vibrating pan»