



HC - SERIES EMMI-100HC



The unit is supplied with lid and basket.



TECHNICAL DATA **EMMI-100HC**

Description:

The Emmi-100HC is the first 10L device in the world where products can be placed along. Thus, the device allows the user to clean a lot of products without contamination.

Dimensions:

External dimensions (L x W x H):
380 x 285 x 295 mm
Inside of Bath (L x W x H): 300 x 245 x 150 mm
Inside of basket (L x W x H): 245 x 180 x 120 mm
Basket mesh size (L x W): 11 x 11 mm

Content:

Tank Capacity: 10 Liters
Working Capacity: 8.5 Liters

Equipment:

Vibrator pan: stainless steel 1.4301
Housing: stainless steel 1.4301
Degree of protection: IP 32
Drain valve: Yes
Timer: 1-60 min / continuous operation
Heating: 20-80 °C
Ultrasonic power controller: 50 / 75 / 100%

Technical data:

Oscillating system:
6PZT wide beam vibratory system
HF Generator: Overload protected, constant
Power operating mode: simple
Half-wave operating frequency: approx. 40 kHz
Ultrasonic input power: 360 watts
Total power in W: 860 Watt
Heating power: 500 Watt
Line voltage: 50 / 60 Hz
Power connection: 220 - 240 V

Accessories:

Stainless steel lid
«noise-resistant»
Stainless Steel Basket
"Protects cleaning material and vibrating pan»