



HC - SERIES EMMI-800HC



The unit is supplied with lid and basket.



TECHNICAL DATA **EMMI-800HC**

Description:

The Emmi-800HC is the largest ultrasonic cleaning device of the Professional Series with drain tap and the first with a Level Meter.

Dimensions:

External dimensions (L x W x H):
770 x 565 x 505 mm
Inside the tub (L x W x H): 600 x 500 x 300 mm
Inside the basket (L x W x H): 560 x 450 x 260 mm
Basket mesh size (L x W): 10 x 10 mm

Capacity:

Bath Capacity: 80 Liters
Working Capacity: 70 Liters

Equipment:

Vibrator pan: stainless steel 1.4301
Housing: stainless steel 1.4301
Degree of protection: IP 20
Drain valve: Yes
Timer switch: 1 - 60 min / continuous operation
Heating: 30-90 °C
Ultrasonic power controller: 50 / 75 / 100%

Technical data:

Oscillating system:
20PZT wide beam vibratory system
HF Generator: Overload protected, constant
Power operating mode: simple
Half-wave operating frequency: approx. 28 kHz
Ultrasonic power: max. 2500 watts
Total wattage: 4500 W
Heating power: 2000 Watt
Line voltage: 50 / 60 Hz
Power connection: 220 - 240 V

Accessories:

Stainless steel lid
«noise-resistant»
Stainless steel basket
«protects cleaning material and vibrating pan»
Basket-inside (L x w x H): 560 x 450 x 260 mm
Basket-mesh size (L x W): 10 x 10 mm