

ADVANTAGES

- + Advanced functionality
- + Modularly expandable
- + Compliant with 21 CFR Part 11

OPTIONS

- LIMS integration
- Suitable with UTS, CIW, DISI, HC6.2, AE
- OPC interface
- Windows authentication via active directory
- WERUM integration

SOFTWARE

PH21 software

For all testing devices and quality assurance

It's time to be 21CFR Part11 compliant. Enjoy high functionality and full data integrity with our software.

The PH21 pharmaceutical quality assurance system enables you to control and evaluate your tablet testers, disintegration testers and weighing machines centrally. Once it is stored in the central product database, you can use product-specific data for all tests on the connected devices.

The PH21 system supports the connection of up to 32 external pharmaceutical testing devices. Such devices include Kraemer Elektronik's well-known UTS tablet testing systems, tablet hardness testers and disintegration testers, as well as weighing systems for in-process control.

For larger applications, the PH21 system can be installed and operated as a client-server application. The entire software is 100% compliant with FDA 21 CFR Part 11 and allows you to use a wide scope of evaluation options for completed tests. Automatic backups in the background guarantee failsafe in-process control.

Numerous interfaces are available for communication and data exchange with external software applications.

Networking example:

