

## STANDARD DEVICE

## DISI-M

Semi-automatic disintegration tester  
with LCD display and individual timer

The DISI-M tablet disintegration tester offers user-friendliness in a tried-and-tested design. Opt for the manual DISI-M disintegration tester as a cost-effective alternative to the automatic DISI-A model. The DISI-M fulfils all EP/USP requirements.

The testing baskets can be attached and removed quickly thanks to the magnetic coupling. A few quick steps is all it takes to dismantle the testing baskets ready for cleaning.

As each station is powered separately, you can test up to four different products at the same time. The softstart function ensures that the testing basket is immersed in a controlled manner.

MODEL	Stations (pieces)	L x W x H (mm)	Weight (kg)
<b>DISI-1M</b>	1	340 x 250 x 750	20
<b>DISI-2M</b>	2	475 x 250 x 750	30
<b>DISI-3M</b>	3	610 x 250 x 750	40
<b>DISI-4M</b>	4	750 x 250 x 750	50



The DISI-M disintegration tester is available with a choice of 1 to 4 measurement stations

## MEASUREMENT PARAMETERS



DISINTEGRATION TIME

## ADVANTAGES

- + Can be upgraded to DISI-A
- + Magnetic coupling of the testing baskets
- + Available with up to four independent measuring stations
- + Stainless steel housing

## OPTIONS

- Testing basket B for tablets from Ø 18 mm
- Testing basket for soft gelatine capsules



Simple & practical: the testing baskets are magnetically coupled



Special testing basket for large tablets